



iTero Scanner Instructions (Two Approaches)

Pre-iTero Scan Instruction:

1. Scan as an iRecord.
2. Do not choose to send to a lab.
3. After you scan your patient's mouth, the scan can be found immediately on MyiTero.com or between approximately fifteen minutes to two days on MyAligntech.com.

Option 1: MyiTero.com Export to STL File Instruction:

1. Log into MyiTero.com.
2. Enter the patients name into the search field (if you can't find it, make sure you are choosing the correct scanner location).
3. Click on the patient's name.
4. Click "Export".
5. An export pop-up box will open and you will see these drop down questions:
 - a. **Export type:** Choose "Solid model with low profile base"
 - b. **Data format:** Choose "Files per arch (arches oriented in occlusion)"
 - c. **File type:** Choose "STL, no color"
 - d. Click "Export".

Note: The STL file is large and may take a minute or two to download. If the export is taking too long, check to make sure that your pop-up blocker is **off** for this site.
6. A window will pop up containing several files associated with your scan. You can copy and paste, drag and drop, or select all, and then save these files to your patient's file.


Option 2: MyAligntech.com Export to STL File Instruction:

1. Log into MyAligntech.com.
2. Click on the ID number for your patient.
3. Click "Export (OrthoCAD 3.5 or higher)".
4. OrthoCAD export pop-up box will open and you will see these drop down questions:
 - a. **Export type:** Choose "Solid model with low profile base"
 - b. **Data format:** Choose "Files per arch (arches oriented in occlusion)"
 - c. **Folder name:** Enter your patient's name in the same format as your chart.
 - d. **File type:** Choose "STL, no Color"
 - e. Click "Export".
5. A dialogue box will pop up asking if you would like to open the exported file's location. Choose "Yes".
6. A window will pop up containing several files associated with your scan. You can copy and paste, drag and drop, or select all, and then save these files to your patient's file.



How to attach files to your patient case once the scans have been exported to your desktop:

1. Login to your DDX Account



2. Select New Case



3. Enter all required patient information and details in the work request form

Portal: SLX Clear Aligner Lab - Cases - Finances - New Task - Search

Schedule a New Case

Patient Information

Provider *	Dr. John Doe	Gender	
First Name		Chart Number *	12345 Test Case
Last Name		Birth date	mm/dd/yyyy

Work Requested

SLX Clear Aligners
Standard, Lite, or Express Package

Remove Procedure

Primary Concern

Chief Complaint *

Case Properties

Impression Type * Treatment Steps * Lab to Recommend *

Arch Development

Archform Width * Optimize Midline * Optimize

A-P Relationship

A-P relationship: Right Canine	Optimize Class I	Overjet / Underbite *	Optimize
A-P relationship: Left Canine	Optimize Class I	Overbite *	Optimize
A-P relationship: Right Molar	Maintain	Crossbite *	Optimize
A-P relationship: Left Molar *	Maintain		



4. Under Case Files, select the upload button to add scan files from your desktop folder

Enclosures

Case Files **Upload** 12 Specified Files Recommended

No files uploaded

Notes

Send Date: 05/06/2021

Emergency

Tags

Coupon Code

Case for SLX Aligner Aligner Lab Uploads

Select Files to Add

DDX File Uploader - DDX - AVG Secure Browser

Open

562021-DDX aligner test

Name	Date modified	Type	Size
562021-DDX aligner test-lowerjaw.stl	5/6/2021 1:36 PM	Meshlayer Docu...	11,195 KB
562021-DDX aligner test-occlusionfirst...	5/6/2021 1:36 PM	Meshlayer Docu...	2,594 KB
562021-DDX aligner test-occlusionsec...	5/6/2021 1:36 PM	Meshlayer Docu...	8,005 KB
562021-DDX aligner test-upperjaw.stl	5/6/2021 1:36 PM	Meshlayer Docu...	13,352 KB

File name: "562021-DDX aligner test-upperjaw.stl" "562021-DDX aligner test-lowerjaw.stl" "562021-DDX aligner test-occlusionsecond.stl" "562021-DDX aligner test-occlusionfirst.stl" All Files (*)

Open Cancel

Select Files to Add

DDX File Uploader - DDX - AVG Secure Browser

ddetstaging.com

Case for SLX Aligner Aligner Lab Uploads

562021-DDX aligner test-lowerjaw.stl	100%	File Uploaded	Patient: Unknown	Specified File: [Dropdown]
562021-DDX aligner test-occlusionfirst.stl	100%	File Uploaded	Patient: Unknown	Specified File: [Dropdown]
562021-DDX aligner test-occlusionsecond.stl	100%	File Uploaded	Patient: Unknown	Specified File: [Dropdown]
562021-DDX aligner test-upperjaw.stl	87.92%	Uploading	Patient: Unknown	Specified File: [Dropdown]

Select Files to Add



5. The uploaded files will appear as thumbnails under the Case Files section. Enter the notes relevant to your patient, check to accept terms & conditions before submitting your aligner case.

Case Files 12 Specified Files Recommended

562021-DDX aligner test-lowerjaw_1.stl
562021-DDX aligner test-occlusionfirst_1.stl
562021-DDX aligner test-occlusionsecond_1.stl
562021-DDX aligner test-upperjaw_1.stl

Notes

Send Date: 05/06/2021

Tags: Try-in, Emergency

Coupon Code

I accept all Terms and Conditions

Submit Case Save Draft

1000 characters remaining

Terms of Use Comments or Suggestions Copyright 2021 Labnet Dental Lab Systems - All Rights Reserved Powered by DDX